


CJACK



AI-Enabled Smart Wearable CPR Device

"Even Saving a Single Life is a Great Victory."





Every Minute Without CPR Costs a Life

Instant Oxygen Loss

Sudden Cardiac Arrest stops the heart — oxygen deprivation to the brain begins within seconds of collapse.

7–10% Drop Per Minute

Survival probability falls 7–10% for every minute CPR is delayed — making speed the difference between life and death.

No Help Nearby

Most cardiac arrests occur outside hospitals, where trained responders are rarely present or accessible in time.

The CPR Gap: Why Bystanders Fail to Act



1. NO KNOWLEDGE

Most bystanders have never received CPR training and don't know where to begin.



REAL-WORLD RESULT

Victims don't receive immediate help. Every minute without CPR reduces survival chances by 7–10%.



7–10%
LESS CHANCE
OF SURVIVAL
EVERY MINUTE



2. FEAR OF ERROR

Hesitation caused by the fear of performing incorrect techniques — or causing harm.



REAL-WORLD RESULT

People freeze or back away, wasting precious time when every second matters.



3. PHYSICAL EXHAUSTION

Manual CPR quality degrades rapidly as fatigue sets in — compressions weaken within minutes.



REAL-WORLD RESULT

Poor quality compressions lead to inadequate blood circulation, reducing the victim's chance of survival.



4. SCALE PROBLEM

One rescuer can only attend one patient at a time — catastrophic in multi-victim scenarios.



REAL-WORLD RESULT

In mass emergencies, many victims don't get help in time, leading to preventable loss of lives.



BRIDGING THE CPR GAP

We need solutions that educate, guide, assist, and scale CPR to empower every bystander to save lives — confidently and effectively.



Existing Solutions Fall Short



MANUAL CPR



REQUIRES TRAINING

Effective CPR requires proper training and practice.



CAUSES FATIGUE

Rescuer fatigue leads to decreased compression rate and depth.



INCONSISTENT COMPRESSION QUALITY

Hard to maintain the correct depth and rate consistently.

CJACK



PORTABLE

Lightweight and easy to carry anywhere.



WEARABLE

Designed like a vest for hands-free operation.



AI-GUIDED

Smart sensors and AI ensure correct depth and rate.



PUBLIC-READY

Built for everyday people and real-world emergencies.



ZERO TRAINING NEEDED

Intuitive to use – anyone can save a life.



LUCAS 3 / AUTOPULSE



BULKY HOSPITAL-GRADE MACHINES

Large and heavy, designed for hospital use.



EXPENSIVE

High cost limits availability for most individuals and organizations.



REQUIRES PROFESSIONAL OPERATION

Needs trained personnel to set up and operate.



NOT PUBLICLY ACCESSIBLE

Not easily available in public places or for bystanders.

The Critical Gap

Hospital-grade mechanical CPR devices like LUCAS 3 and AutoPulse deliver consistent compressions — but they are designed for clinical environments, not public spaces.

Manual CPR demands training and tires rescuers within minutes. Neither solution addresses the everyday emergency in a mall, office, or street.

✗ No portable, public-ready, automated CPR solution currently exists for everyday emergencies.



THE FUTURE OF CPR IS ACCESSIBLE.

CJACK bridges the gap between manual CPR and expensive machines, making high-quality CPR possible for everyone, everywhere.



Introducing CJACK

CJACK is an **AI-enabled, fully wearable CPR device** designed for immediate deployment in public spaces — delivering standardized, uninterrupted chest compressions until emergency medical professionals arrive.

Wearable Jacket

Fits any adult patient — deployed in seconds from a wall-mounted station.

AI Voice Assistant

Guides any bystander step-by-step — zero prior CPR training required.

Auto Emergency Call

Automatically dials **108** emergency services the moment the device is activated.

WHEN IS CJACK REQUIRED?



Cardiac
Arrest



Respiratory
Arrest



Drowning



Electric
Shock



Drug
Overdose



Severe
Trauma



Airway
Obstruction

RAPID DEPLOYMENT – 30 SECONDS

FAST. SIMPLE. LIFE-SAVING.

1 OPEN THE JACKET



 Open and prepare the CJACK for deployment.

2 PLACE JACKET ON PATIENT



 Place the jacket over the patient and secure the shoulder straps.

3 ALIGN AT STERNUM



 Align the center of the device precisely over the lower sternum.

4 LOCK THE STRAPS



 Tighten and lock all straps for a secure and stable fit.

5 PRESS START BUTTON



 Press the START button to activate the system.

6 CPR & OXYGEN ASSIST STARTS AUTOMATICALLY



 Automated compressions and oxygen assist begin immediately.



REAL-TIME
FEEDBACK



ADAPTIVE
COMPRESSION



DEPTH
TRACKING



WIRELESS
MONITORING



GPS
TRACKING



MULTI-USER
SYNC



IP65
WATERPROOF



LONG BATTERY
LIFE



BLUETOOTH



WIFI



4G
LTE

Design of CJACK



Structure of CJACK:



DEPLOY ANYWHERE – LIKE AN AED



AMBULANCE



SCHOOL



GYM



TEMPLE / RELIGIOUS CENTER



COLLEGE / HOSTEL



FACTORY



HOME



RAILWAY STATION



BUS STAND



SPORTS EVENTS



RURAL AREAS



PUBLIC PLACES

MULTIPLE VICTIM SCENARIO

Faster Help. More Survivors. Better Outcomes.

✗ WITHOUT CJACK

5 victims collapse. Only 2 rescuers available.



2 patients receive CPR



3 patients receive **no help at all**



Rescuers are physically exhausted and unable to respond further

✓ WITH CJACK

Rescuers deploy jackets **on all 5 patients** simultaneously.



All 5 patients receive **continuous, automated CPR**



Rescuers are freed to coordinate, communicate, and assist further



Emergency services receive an automatic alert the moment each device activates



CJACK MULTIPLIES IMPACT.

One rescuer can save more lives with CJACK.



2 HELPED



5 HELPED

■ AHA-Compliant CPR Standards



Built to AHA Guidelines

CJACK's compression engine is engineered to meet the precise parameters set by the **American Heart Association (AHA)** — the gold standard in resuscitation science.

© 2020 CJACK Inc. All rights reserved.

1

100–120 / min

Compression Rate

© 2020 CJACK Inc. All rights reserved.

2

5–6 cm

Compression Depth

© 2020 CJACK Inc. All rights reserved.

3

Lower Sternum

Center of chest — precise positioning every time

What Makes CJACK Unique



Wearable & Portable

Unlike hospital machines, CJACK can be installed **anywhere** — malls, airports, offices, stadiums — mounted on walls like a fire extinguisher.



Public Emergency Device

Mounted with a clear instruction board, CJACK is as familiar and accessible as a fire extinguisher — a true public safety tool.




AI-Guided Operation

Simple voice instructions walk any bystander through the process. **Zero prior training required** — designed for the untrained public.



Multiple Patient Support

Each jacket operates **independently** — rescuers can simultaneously deploy devices across multiple victims during mass casualty events.



Selling price and its types

 SALES TYPE	 PRICING
 Individual Purchase	₹29,999
 Bulk Hospital Orders	₹27,000 per unit
 Industrial Safety Orders	₹26,500 per unit
 Government/NGO Partnership	Negotiable pricing

 PRODUCT VERSION	 SELLING PRICE (₹)
 BASIC Basic Model	₹21,000 – ₹23,000
 STD Standard Smart Model	₹25,000 – ₹25,000
 AI Advanced AI Emergency Model	₹18,000 – ₹35,000

Key Features at a Glance



AI Voice Guidance



Fully Automated CPR



Wearable & Portable



Auto 108 Calling



Emergency Stop




Fatigue-Free Operation



Public Wall Installation



AHA-Standard Compressions



AI-ENABLED

AHA-COMPLIANT

PUBLIC SAFETY INNOVATION

CJACK

Because Every Second Counts.

From the moment someone collapses to the moment help arrives — CJACK bridges the gap that costs thousands of lives every year. Wearable. Automated. AI-Guided. Ready for anyone, anywhere.